




# One-Stop Solution for Western Blotting

Separate



## SurePAGE™ Pre-cast Gels

Fast run times of 20+ Min  
High performance  
Free running buffer



## Protein Standards

Comprehensive selection & Broad  
Molecular Weight Range

Transfer  
&  
Stain



## eBlot™ Protein Transfer Systems

Efficient wet transfer of mini and  
midi gels in 5-15 minutes



## eStain® L1 Protein Staining System

Coomassie Blue or Ponceau Staining  
in 10 Min

Wash  
&  
Incubate



## eZwest™ Lite Automated Western Device

Automated Western Blotting device completing the entire  
process of blocking, antibody incubation and washing steps.

Image  
&  
Quantify



## THE™ Elite Antibodies

Precise and effective antibodies that meet the most  
demanding research requirements.

# Ordering Information

## SurePAGE Pre-Cast Gels

Cat No.	Product Name	Quantity
M00661	SurePAGE™, Bis-Tris, 10×8, 8%, 10 wells	10 PK
M00662	SurePAGE™, Bis-Tris, 10×8, 8%, 12 wells	10 PK
M00663	SurePAGE™, Bis-Tris, 10×8, 8%, 15 wells	10 PK
M00664	SurePAGE™, Bis-Tris, 10×8, 10%, 10 wells	10 PK
M00665	SurePAGE™, Bis-Tris, 10×8, 10%, 12 wells	10 PK
M00666	SurePAGE™, Bis-Tris, 10×8, 10%, 15 wells	10 PK
M00667	SurePAGE™, Bis-Tris, 10×8, 12%, 10 wells	10 PK
M00668	SurePAGE™, Bis-Tris, 10×8, 12%, 12 wells	10 PK
M00669	SurePAGE™, Bis-Tris, 10×8, 12%, 15 wells	10 PK
M00652	SurePAGE™, Bis-Tris, 10×8, 4-12%, 10 wells	10 PK
M00653	SurePAGE™, Bis-Tris, 10×8, 4-12%, 12 wells	10 PK
M00654	SurePAGE™, Bis-Tris, 10×8, 4-12%, 15 wells	10 PK
M00655	SurePAGE™, Bis-Tris, 10×8, 4-20%, 10 wells	10 PK
M00656	SurePAGE™, Bis-Tris, 10×8, 4-20%, 12 wells	10 PK
M00657	SurePAGE™, Bis-Tris, 10×8, 4-20%, 15 wells	10 PK
M00658	SurePAGE™, Bis-Tris, 10×8, 8-16%, 10 wells	10 PK
M00659	SurePAGE™, Bis-Tris, 10×8, 8-16%, 12 wells	10 PK
M00660	SurePAGE™, Bis-Tris, 10×8, 8-16%, 15 wells	10 PK

## Protein Standards

Cat No.	Product Name	Quantity
M00624	Broad Multi Color Pre-Stained Protein Standard	250 uL, 1250 uL
M00516	PAGE-MASTER Protein Standard	500 µL
MM1397-500	PAGE-MASTER Protein Standard Plus	500 µL
M00521	WB-Master Protein Standard	250 µL

## eBlot Fast Transfer System

Cat No.	Product Name
L00686	eBlot L1 Fast Transfer Device
L00724	NC eBlot L1 Transfer Kit (Excludes NC Membranes)
L00726	PVDF eBlot L1 Transfer Kit (Excludes PVDF Membranes)
L00727	Complete eBlot L1 Transfer Kit (PVDF)
L00981	eBlot™ L2 Fast Protein Transfer Device
L01015-30	eBlot L2 Transfer Basic Kit, midi
L01015-30 (1) + L01067 (1)	eBlot L2 Midi Package with PVDF Midi Membranes
L01015-30 (1) + L01066 (2)	eBlot L2 Midi Package with PVDF Mini Membranes

## eStain L1 Protein Staining System

Cat No.	Product Name
L00657	eStain L1 Protein Staining Device
L00753	eStain L1C Protein Staining Kit

## eZwest Automated Western Device

Cat No.	Product Name
L00816	eZwest Automated Western Device